



COMPARISON OF BODYSTAT

	Bodystat®1500 Single Frequency Adults	Bodystat®1500MDD Dual Frequency Adults & Children
DISPLAY OPTIONS		
Fat Mass (FM)* [% & kg]	✓	✓
Lean Mass (FFM)* [% & kg]*	✓	✓
Dry Lean Mass* [kg]	✓	✓
Total Weight (max 300 kg)	✓	✓
Total Body Water (TBW) * [% & lt]	✓	✓
Basal Metabolic Rate (BMR)* [kcal] (Based on FFM vs Total Body Weight)	✓	✓
BMR/Body Weight* [kcal/kg]	✓	✓
Activity Energy Requirement*	✓	✓
Waist/Hip Ratio (WHR)	✓	✓
Body Mass Index (BMI)	✓	✓
Body Fat Mass Index (BFMI)*		✓
Fat Free Mass Index (FFMI)*		✓
Wellness Marker™		✓
Extra-cellular Water (ECW)*		
Intra-cellular Water (ICW)*		
Nutrition Index (ECW/TBW)*		
Body Cell Mass (BCM)*		
3rd Space Water*		






COMPARISON OF BODYSTAT

Illness Marker™		
<u>Dehydration Index</u>		
<u>Segmental Impedance Analysis</u>		
Frequency Measured	50 kHz	5/50 kHz
Impedance (Range 20 - 1,000 ohms)	✓	✓
Resistance (50 kHz)		✓
Reactance (50 kHz)		✓
<u>Phase Angle</u> (50 kHz)		✓
Normal Ranges **	✓	✓

ADDITIONAL FEATURES

Automatic 100 Test Memory Storage for subsequent recall	✓	✓
Batteries Duracell AA-size 1.5V	6	6
Real Time Clock records exact time of measurement		✓
Bluetooth Communications (Software/Printer)	Only supplied with Software/Printer	✓
European Foreign Languages	✓	✓
Portable Thermal Bodystat Printer	Optional	Optional
Warranty 12 months (excludes lead wires)	✓	✓
Lead Wires Shielded (preventing electrical interference)	✓	✓
Enclosures Shielded (preventing electrical interference)		
Electrodes Used	Short	Long



COMPARISON OF BODYSTAT

AGE RANGE FOR BODY COMPOSITION

6 - 17 years	Available in Wellness and Body Manager Software	✓
18 - 70 years	✓	✓
> 70 years	Available in Body Manager Software	Included in Body Manager Software

BODYSTAT WINDOWS SOFTWARE

Wellness - up to 17 page educational report	Optional	Optional
Body Manager - includes Trend Reports	Optional	Included
QuadScan - includes Hydration Analysis		

